



SPORT RIDER

drive

A MILESTONE IN
SCOOTER DESIGN
BUILT TO PUSH
THE LIMITS





THE STYLISH NEW SPORT RIDER WAS
INSPIRED BY SOME OF THE LEADING
DESIGNS IN THE MOTORCYCLE
INDUSTRY. THIS MOVES AWAY
FROM THE TRADITIONAL TYPE
OF SCOOTER, CREATING
DESIRABILITY AND EXCITEMENT
TO THE DRIVING EXPERIENCE.






THE SPORT RIDER IS A MODERN DESIGN 8MPH SCOOTER THAT INCLUDES ALL THE TRIED AND TESTED BENEFITS OF THE DRIVE RANGE. THIS EASY TO OPERATE VEHICLE PROVIDES STATE OF THE ART PERFORMANCE AND WILL GIVE YEARS OF DEPENDABLE SERVICE. USING TECHNOLOGY BORROWED FROM THE AUTOMOTIVE INDUSTRY THE ALL ROUND SUSPENSION PROVIDES EXCELLENT RIDE QUALITY AND HANDLING, MAKING A CRUISE ON THE SPORT RIDER BOTH EXHILARATING AND COMFORTABLE. THE SPORT RIDER CARRIES THE DRIVE BRAND, THIS ENSURES A FOCUS ON COMFORT, STYLE, PERFORMANCE AND RELIABILITY.



THE SPORT RIDER IN DETAIL

ALL ROUND LIGHTING INCLUDING A LARGE HEADLIGHT, FRONT & REAR INDICATORS AND REAR BRAKE LIGHTS



TELESCOPIC FRONT
MOTORCYCLE
SUSPENSION
FOR OPTIMAL
HANDLING



MOTORCYCLE STYLE WHEELS WITH CHROME MUD GUARDS



TWIN DIAL DIGITAL DASH & TRIP COMPUTER WITH ODOMETER

SPORT RIDER FEATURES



EASY TO USE CONTROLS,
TWIN DIAL DIGITAL
DASH INCLUDING A TRIP
COMPUTER WITH ODOMETER.
FULLY AUTOMATIC EASY TO
USE "SELF DIAGNOSTIC"
FAULT-FINDING SYSTEM.

ALL ROUND LIGHTING INCLUDING
A LARGE HEADLIGHT, FRONT AND
REAR INDICATORS AND REAR
BRAKE LIGHTS. TELESCOPIC FRONT
MOTORCYCLE SUSPENSION
FOR OPTIMAL HANDLING.
EXCELLENT MANOEUVRABILITY
AND TURNING CIRCLE.

COMFORTABLE ROTATING
SPORT STYLE SEAT WITH
SLIDING AND RECLINING
FACILITY. WIDTH ADJUSTABLE
ARMRESTS, HEADREST AND
A LARGE SEAT BACK
STORAGE POCKET.

REAR BUMPER PROVIDES
ADDITIONAL PROTECTION.
4MPH PAVEMENT AND 8MPH
MAXIMUM ROAD SPEED SETTINGS
WITH A 31 MILE RANGE*. LARGE
ACCOMMODATING FOOTWELL
WITH ANTI SLIP RUBBER MAT.

*SUBJECT TO CONDITIONS OF USE


MOTORCYCLE STYLE
WHEELS WITH
CHROME MUD GUARDS
AND ANTI TIP
WHEELS.

THROTTLE LEVER WITH
AUTOMATIC SAFE-BRAKING
SYSTEM PLUS AN ADDITIONAL
HAND BRAKE. ANGLE
ADJUSTABLE HANDLE
BARS AND KEY
IGNITION FOR SECURITY.

SPECIFICATION & CONTACT

LENGTH	165CM (65")
WIDTH	78CM (31")
HEIGHT	112CM (44")
TURNING RADIUS	157CM (62")
GROUND CLEARANCE	12CM (5")
MAX RANGE	31 MILES
WEIGHT LIMIT	180KG (28ST)
MAX SPEED	8MPH
BATTERIES	100AH 12V X 2
BRAKE SYSTEM	ELECTROMAGNETIC
CLIMBING ANGLE	12°
CHARGING TIME	12 HRS
CONTROLLER	160 AMP/DYNAMIC
DISPLAY	DIGITAL
COLOUR	BLACK & CHROME
LEG ROOM	51CM (20")
SEAT TYPE	SPORTS STYLE
SUSPENSION	FULL FRONT & REAR
PANNIER	OPTIONAL EXTRA
MIRRORS	LEFT & RIGHT AS STANDARD

FOR MORE INFORMATION PLEASE CONTACT YOUR LOCAL DEALER:



COMFORT
STYLE
PERFORMANCE
RELIABILITY

