

CONTACT YOUR LOCAL DEALER FOR MORE INFORMATION



Royale
SCOOTERS

drive[™]
WWW.ROYALESCOOTERS.CO.UK

drive[™]

ROYALE

3
&
4

THE RESULT

EXCEPTIONAL

8MPH*

DRIVING

EXPERIENCE

ERGONOMIC SPORT SEAT

FULLY ACTIVE
INDEPENDENT
FRONT & REAR
SUSPENSION

EXTRA-LARGE
16" REAR WHEELS
PNEUMATIC TYRES

ROYALESCOOTERS.CO.UK

ROYALE

LCD DISPLAY

ALLOY
WHEELS

3
&
4

STYLISH
DESIGN

INTELLIGENT DIAGNOSTICS

MAX RANGE
32 MILES

SPEEDO


ODOMETER

POWER
INDICATOR

70AMP BATTERIES



**THE DRIVE ROYALE
3 & 4 ARE TWO
OF THE FINEST
MOBILITY SCOOTERS
AVAILABLE TODAY.**



FEATURES

- MAX RANGE OF 32 MILE* / TOP SPEED OF 80PH*
- 18KG (28ST) WEIGHT LIMIT / ALL ROUND SUSPENSION / DIGITAL DASH / SLOING AND ROTATING SEAT / HIGH LEVEL INDICATORS / STYLISH DESIGN / 70AH BATTERIES / 4 COLOURS AVAILABLE
- AUTOMATIC SPEED REDUCTION WHEN STEERING
- LARGE 16" DIAMETER WHEELS / ARMRESTS ARE ANGLE ADJUSTABLE / TILLER ADJUSTMENT BY LEVER
- TWIN REAR VIEW MIRRORS / ROAD & PAVEMENT LEGAL / COMFORTABLE SEAT WITH RECLINING FACILITY / HEIGHT ADJUSTABLE HEADREST / OPTIONAL REAR PANNIER

BOTH MODELS GIVE AN IMPRESSIVE PERFORMANCE;

POWERED BY 70AH BATTERIES THEY HAVE AN EXCELLENT MAXIMUM RANGE OF 32 MILES. THESE STYLISH SCOOTERS ALSO FEATURE A MODERN ILLUMINATED LCD DISPLAY, INDICATING SPEED, POWER, CLOCK, ODOMETER, TEMPERATURE, REVERSE LIGHT AND INTELLIGENT DIAGNOSTICS. WITH FULLY ACTIVE INDEPENDENT FRONT AND REAR SUSPENSION, EXTRA-LARGE 16" REAR WHEELS, PNEUMATIC TYRES AND AN ERGONOMIC SPORT SEAT, THE RESULT IS SIMPLE, AN EXCEPTIONAL DRIVING EXPERIENCE.

COLOUR OPTIONS

RED	WHITE
BLACK	SILVER

ROYALE 4

SPECIFICATION

CODE

HW002

LENGTH

158CM (62")

HEIGHT

135CM (53")

WIDTH

73CM (29")

MAX CAPACITY

185KG (28ST)

WEIGHT

149KG (330LBS)

MAXIMUM GRADIENT

10°

GROUND CLEARANCE

150MM (6")



ROYALE 3

SPECIFICATION

CODE

HW001

LENGTH

158CM (62")

HEIGHT

135CM (53")

WIDTH

73CM (29")

MAX CAPACITY

185KG (29ST)

WEIGHT

149KG (330LBS)

MAXIMUM GRADIENT

10°

GROUND CLEARANCE

150MM (6")